

Batterie standard européenne

(mise à jour 2023)

	Allergènes	Concentration
1	BICHROMATE DE POTASSIUM	0.5% vas
2	P-PHENYLENEDIAMINE	1% vas
3	THIURAM MIX	1% vas
4	NEOMYCINE SULFATE	20% vas
5	CHLORURE DE COBALT	1% vas
6	FORMALDEHYDE	2% aq ¹
7	COLOPHONIUM	20% vas
8	2-HYDROXYETHYL METHACRYLATE (HEMA)	2% vas ¹
9	MYROXYLON PEREIRAE	25% vas
10	N-ISOPROPYL-N-PHENYL-4-PHENYLENEDIAMINE (IPPD)	0.1% vas
11	LANOLIN ALCOHOL	30% vas
12	MERCAPTO MIX	2% vas
13	EPOXY RESIN	1% vas
14	PARABEN MIX	16% vas
15	4- TERT-BUTYL PHENOL FORMALDEHYDE RESIN	1% vas
16	FRAGRANCE MIX 1	8% vas
17	SESQUITERPENE LACTONE MIX	0.1% vas
18	METABISULFITE DE SODIUM	1% vas
19	NICKEL SULFATE HEXAHYDRATE	5 % vas
20	METHYLISOTHIAZOLINONE + METHYLCHLOROISOTHIAZOLINONE	0.02% aq ¹
21	2-MERCAPTOBENZOTHAZOLE	2% vas
22	METHYLDIBROMO GLUTARONITRILE	0.5% vas
23	BUDESONIDE	0.01% vas
24	TIXOCORTOL-21-PIVALATE	0.1% vas
25	TEXTILE DYE MIX	6.6 % vas
26	FRAGRANCE MIX II	14% vas
27	HYDROXYISOHEXYL 3-CYCLOHEXENE CARBOXALDEHYDE	5% vas
28	PROPOLIS	10% vas
29	METHYLISOTHIAZOLINONE	0.2% aq ¹
30	CAINE MIX 3	10% vas
31	DECYL GLUCOSIDE	5 % vas
32	BENZISOTHIAZOLINONE	0..1% vas

¹ Remplir la chambre au moment de l'application